# Hera

## HERAVEST® Onyx

DE	Gebrauchsanweisung .....	3
GB	Instructions for use .....	6
FR	Mode d'emploi .....	9
ES	Instrucciones de uso .....	12
IT	Istruzioni per l'uso .....	15
NL	Gebruiksaanwijzing .....	18
SE	Bruksanvisning .....	21

**KULZER**  
MITSUI CHEMICALS GROUP

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**Graphite-free, phosphate-bonded, precision investment material for “speed” and program-med preheating crown and bridgework, especially for non-precious alloys**

### Storage

Keep the powder in a cool, dry place

Do not store the mixing liquid at less than 5°C as it is sensitive to frost.

### Working temperature

approx. 22°C (Room temperature)

### Mixing ratio

100 g Powder – 22 ml Liquid

60 g Powder – 13 ml Liquid

### Mixing ratio for 160 g Sachet:

Powder 160 g, Liquid 35 ml

For accurate dispensing please use our dispensing syringes.

Liquid concentration (examples)		
Liquid concentration for 1000 ml liquid		
Concentration	Concentrate (ml)	Distilled water (ml)
50 %	500	500
65 %	650	350
75 %	750	250
95 %	950	50
100 %	1000	–

Concentrations						
Concentration %	Sachet 1 x 160 g		Sachet 2 x 160 g (320 g)		Sachet 3 x 160 g (480 g)	
	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)
100	35	–	70	–	105	–
95	33	2	66	4	99	6
90	32	3	64	6	96	9
85	30	5	60	10	90	15
80	28	7	56	14	84	21
75	26	9	52	18	78	27
70	25	10	50	20	75	30
65	23	12	46	24	69	36
60	21	14	42	28	63	42
55	19	16	38	32	57	48
50	17,5	17,5	35	35	52,5	52,5

► Expansion control (%)				
Indication/ Pattern material		alloy	Bonding alloys Nickel chrome molybdenum (NiCrMo)	Bonding alloys Cobalt chrome molybdenum (CoCrMo)
		1. Crown and bridgework	- Wax	
2. Outer copings and small-span bridges	- Pattern Resin		65 – 75	75 – 85
3. Outer telescopes	- Pattern Resin		–	90 – 100
4. Investment Special BS Liquid 3 can be used to obtain greater expansion. Mixing ratio: 100 g Powder - 22 ml Liquid				

### Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

### ► Mixing in vacuum

60 seconds

Place in vacuum for 15 seconds but without switching the paddles on. Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

### Working time

Approx. 7 min. at ap-prox. 22°C (Room temperature)

The working time depends on the room temperature. Warmth curtails the working time.

► Suitable ring systems	
Use <b>only steel rings</b> with cellulose ring liners	<b>Casting ring – X1 and X3 steel:</b> We recommend using <b>one</b> dry cellulose ring liner with Vaseline. <b>Casting ring – X6 and X9 steel:</b> We recommend using <b>two</b> dry cellulose ring liners with Vaseline.

### Investing

Fill the ring with investment: only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.

### Placing the rings in a preheated furnace

20 minutes after start of mix. Roughen the top of the mould before placing it in the furnace!

### Preheating

► Speed preheating					
Set-up temp. (°C)	Final temp. (°C)	Hold-times (min.) (commences once the final temp. is attained again)			
		X1	X3	X6	X9
850-950	max. 950	30	30	45	60

\*The holding time should be extended if several rings are to be heated at the same time.

▶ Programmed preheating						
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold-times (min.)			
			X1	X3	X6	X9
1	approx. 23°C room temp. - 270°C	8	30	30	45	60
2	270°C - 580°C	8	30	30	45	60
3	580°C - final temp. (950°C)	8	30	30	45	60

### ⚠ Safety Instructions

- Do not open the furnace during preheating as this would cause aggressive combustion.
- The preheat furnace should only be filled to half of its maximum capacity.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

### ⚠ Changes

Revision marking April 2008 / Changes in relation to issue 02/2005.

Expansion control, vacuum stirring time, suitable casting ring systems, preheating process

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

The figures given by our research and development department are guidelines and can be individually adjusted.

www.kulzer.com

We appreciate your comments, feedback, and suggestions.

**Contact in the countries** and more information are **available by the displayed QR code** or at our website **www.kulzer.com**

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Kulzer GmbH