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Made in Germany

KULZER
MITSUBI CHEMICALS GROUP

Preciano® conductive silver varnish

Hera

Hera Preciano conductive silver varnish is a thin consistency, rapid drying lacquer containing silver. It is used for producing electro-conductive surfaces for electroforming restorations. It does not react with Hera Preciano Bio Gold solution universal.

Application: Shake conductive silver varnish carefully before use until the mixing balls can be heard. Place a small amount in the dropper lid and use a Hera Preciano special brush to apply a thin uniform, masking coat of varnish without allowing it to puddle. The surface being electroformed must have been sealed in advance with plaster hardener and allowed to dry – the contact rod, which is adhered into place, must also be coated. All inner telescopes must be kept grease-free and dust-free before applying hardener. Only use minute amounts of varnish to touch up defective areas in the coating. Allow to dry for at least 30 minutes at room temperature. After removing conductive silver varnish from the bottle, cover it immediately with foil (Adapta) etc. to prevent the varnish drying out. Once finished, close the bottle tightly and clean the Hera Preciano special brush with brush cleaner.

Once electroforming has been completed, the conductive silver varnish can be dissolved in an ultrasonic cleaner (max. temperature: 40°C (104°F)) using Hera Preciano conductive silver varnish remover.

Please note: As thick conductive silver varnish may lead to uneven coatings, discontinue using it. Please do not thin!

For further information on the use of conductive silver lacquer for telescopic crowns and bar-borne restorations, please refer to currently valid Hera Preciano instructions for use.

Storage: Seal the container airtight immediately after use, place it in the packaging and keep it in a well ventilated place, protected from sunlight.

Shelf-life: See imprint on paper box and bottle.

Dated: 2017-07

Hazard-determining components: 1-methoxypropan-2-ol.

Hazard statements

Highly flammable liquid and vapour. May cause drowsiness or dizziness.

Precautionary statements

Keep away from heat / sparks / open flames / hot surfaces. – No smoking. Avoid breathing dust / fume / gas / mist / vapours / spray.

Danger

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor / physician if you feel unwell. Dispose of contents / container in accordance with local / regional / national / international regulations.

Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)

Made in Germany

KULZER
MITSUI CHEMICALS GROUP